

### Formulas used:

The reference distance is calculated as follows:

$$r_{ref} = \frac{res_x^0 * res_y^0 * res_z^0}{\sqrt{(res_y^0 * res_z^0 * \cos\Phi * \sin\theta)^2 + (res_x^0 * res_z^0 * \sin\Phi * \sin\theta)^2 + (res_x^0 * res_y^0 * \cos\theta)^2}}$$

$$\Phi = \arccos \frac{x_B - x_A}{\sqrt{(x_B - x_A)^2 + (y_B - y_A)^2}}$$

$$\theta = \arccos \frac{z_B - z_A}{\sqrt{(x_B - x_A)^2 + (y_B - y_A)^2 + (z_B - z_A)^2}}$$

$$NA \leq 0.7 \quad res_{x,y}^0 = \frac{0.377 * \lambda_{ex}}{NA} \quad res_z^0 = \frac{0.626 * \lambda_{ex}}{n - \sqrt{n^2 - NA^2}}$$

$$NA > 0.7 \quad res_{x,y}^0 = \frac{0.383 * \lambda_{ex}}{NA^{0.91}} \quad res_z^0 = \frac{0.626 * \lambda_{ex}}{n - \sqrt{n^2 - NA^2}}$$

$x_A, y_A, z_A$  and  $x_B, y_B, z_B$  are the bead coordinates in channel A and B respectively, NA: numerical aperture,  $\lambda_{ex}$ : excitation wavelength, n: refractive index of the lens immersion & mounting media.

Lateral ( $res_{x,y}^0$ ) and axial ( $res_z^0$ ) theoretical resolution values used for single point scanning multiphoton excitation microscopes are calculated as defined in Zipfel, W.R., Williams, R.M. & Webb, W.W. Nonlinear magic: multiphoton microscopy in the biosciences, Nature Biotechnology 21, 1369–1377 (2003).

Compliance with the Shannon-Nyquist criterion for k photon excitation uses the following formulas for Shannon-Nyquist distances calculation (and a fixed k value of 2):

$$\alpha = \arcsin\left(\frac{NA}{n}\right)$$

$$\Delta_{x,y} = \frac{\lambda_{ex}}{4.k.NA} \quad \Delta_z = \frac{\lambda_{ex}}{4.n.(1 - \cos(\alpha))}$$